# Lyphochek® Coagulation Control

## Levels 1, 2 and 3

<b>REF</b>	<b>744</b>	Level 1	12 x 1 mL	<b>CE</b>	<b>IVD</b>	<b>EXP</b> 2013-06-30	<b>LOT</b> 78200	Level 1	78201
	<b>745</b>	Level 2	12 x 1 mL					Level 2	78202
	<b>746</b>	Level 3	12 x 1 mL					Level 3	78203
	<b>745X</b>	Trilevel MiniPak	3 x 1 mL						

### ENGLISH

#### INTENDED USE

Lyphochek Coagulation Control is intended for use as a quality control plasma to monitor the precision of citrated coagulation systems.

#### SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. Three levels of control are available to allow performance monitoring within the clinical range.

#### REAGENT

This product is prepared from human plasma, with added purified biochemicals and preservatives. The control is provided in lyophilized form for increased stability.

#### STORAGE AND STABILITY

Using a volumetric pipet, reconstitute each vial with 1.0 mL of distilled or deionized water. Replace the stopper and allow the control to stand for 10 to 15 minutes, swirling occasionally. Before sampling, gently invert the vial several times to ensure homogeneity.

capped, all analytes will be stable at:

- 2 to 25°C for 48 hours
- or -20 to -70°C for 30 days
- or -20 to -70°C with three freeze-thaw cycles within 7 days.

This product is shipped under ambient conditions.

#### RECONSTITUTION

Using a volumetric pipet, reconstitute each vial with 1.0 mL of distilled or deionized water. Replace the stopper and allow the control to stand for 10 to 15 minutes, swirling occasionally. Before sampling, gently invert the vial several times to ensure homogeneity.

#### PROCEDURE

This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit, or reagent being used.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.

#### LIMITATIONS

- This product should not be used past the expiration date.
- If there is evidence of microbial contamination or excessive turbidity in the reconstituted product, discard the vial.
- This product is not intended for use as a standard.

#### ASSIGNMENT OF VALUES

The mean values printed in this insert were derived from replicate analyses and are specific for this lot of product. The tests listed were performed by the manufacturer and/or independent laboratories using manufacturer supported reagents and a representative sampling of this lot of control. Individual laboratory means should fall within the corresponding acceptable range; however, laboratory means may vary from the listed values during the life of this control. Variations over time and between laboratories may be caused by differences in laboratory technique, instrumentation and reagents, or by manufacturer test method modifications. It is recommended that each laboratory establish its own means and acceptable ranges and use those provided only as guides.

Refer to [www.qcnet.com](http://www.qcnet.com) for insert update information.

#### SPECIFIC PERFORMANCE CHARACTERISTICS

This product is a freeze-dried product manufactured under rigid quality control standards. To obtain consistent vial-to-vial assay values, the control requires proper storage and handling as described.

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**ENGLISH**

**Biological source material. Treat as potentially infectious.**

Each human donor unit used to manufacture this control was tested by FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), antibody to Hepatitis C (HCV) and antibody to HIV-1/HIV-2. This product may also contain other human source materials for which there are no approved tests. In accordance with good laboratory practice, all human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

Material Safety Data Sheet (MSDS) available for professional users on [www.bio-rad.com](http://www.bio-rad.com).

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ENGLISH

- (1) Representative International Normalized Ratio (INR) values determined using thromboplastin with ISI of 0.75 – 2.25.
- (2) Low Fibrinogen levels may affect clotting times on various instruments. Please refer to instrument manual.
- ▲ Data not available at the time of printing. Please inquire.
- § The data required to establish the means and acceptable ranges for this assay were not obtained due to limited assignment participation. If your facility is interested in participating in the Value Assignment Program for this assay, please contact your local Bio-Rad Sales or Technical Services Group.
- ❖ INTERNATIONAL USE ONLY - The following section contains data for methods that are not available for diagnostic use in the United States.

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## **- INTERNATIONAL USE ONLY -**

**The following section contains data for methods that are not available for diagnostic use in the United States. ❖**

INTERNATIONAL USE ONLY -  
The following section contains data for methods that are not available for diagnostic use in the United States. ❖

METHOD							
Methode // Méthode // Metodo // Método // Método // Metod // Metode							
	Units	Level 1 - 78201		Level 2 - 78202		Level 3 - 78203	
		Mean	Range	Mean	Range	Mean	Range
PROTHROMBIN TIME (PT)							
Axis Shield Nycotest PT	Seconds	\$		\$		\$	
Diagnostica Stago / Roche STA Hepato Quick	Seconds	26.1	20.9 – 31.4	65.9	52.7 – 79.0	81.3	65.1 – 97.6
Diagon Dia-PT	Seconds	▲		▲		▲	
Fortress Diagnostics PT-HS	Seconds	▲		▲		▲	
Axis Shield Nycotest PT	INR (1)	\$		\$		\$	
Diagnostica Stago / Roche STA Hepato Quick	INR (1)	0.960	0.768 – 1.15	2.10	1.68 – 2.52	2.54	2.03 – 3.05
Diagon Dia-PT	INR (1)	▲		▲		▲	
Fortress Diagnostics PT-HS	INR (1)	▲		▲		▲	
Axis Shield Nycotest PT (Calculated)	%	\$		\$		\$	
Diagnostica Stago / Roche STA Hepato Quick	%	110	88.1 – 132	27.8	22.2 – 33.4	21.4	17.1 – 25.7
Diagon Dia-PT	%	▲		▲		▲	
Fortress Diagnostics PT-HS	%	▲		▲		▲	
ACTIVATED PARTIAL THROMBOPLASTIN TIME (APTT)							
Fortress Diagnostics APTT	Seconds	▲		▲		▲	
FIBRINOGEN (2)							
Fortress Diagnostics Fibrinogen (Quantitative)	mg/dL	▲		▲		▲	
	g/L	▲		▲		▲	

## NOTES

FOR REFERENCE USE ONLY

#### Autoimmune

Liquichek™ ANA Controls  
Centromere Pattern  
Homogeneous Pattern  
Mitotic Spindle Pattern  
Nucleolar Pattern  
Speckled Pattern  
Liquichek™ ANA Controls Set, Positive (Homogeneous, Speckled, Centromere & Nucleolar Patterns)  
Liquichek™ Anti-Mitochondrial Control  
Liquichek™ Anti-nDNA Control  
Liquichek™ Anti-RNP Control  
Liquichek™ Anti-Scl-70 Control  
Liquichek™ Anti-Sm Control  
Liquichek™ Anti-Smooth Muscle Control  
Liquichek™ Anti-SS-A Control  
Liquichek™ Anti-SS-B Control  
Liquichek™ Autoimmune Negative Control

#### Blood Gas

Liquichek™ Blood Gas Control  
Liquichek™ Blood Gas Plus CO-Oximeter Control (Bayer 800 Series)  
Liquichek™ Blood Gas Plus E Control  
Liquichek™ Blood Gas Plus EGL Control

#### Cardiac Assessment

Liquichek™ Cardiac Markers Control LT  
Liquichek™ Cardiac Markers Control LT, Level Low  
Liquichek™ Cardiac Markers Plus Control  
Liquichek™ Cardiac Markers Plus Control LT  
Liquichek™ CK/LD Isoenzyme Control  
Liquichek™ Homocysteine Control

#### Chemistry

Lyphocheck® Assayed Chemistry Control  
Liquid Assayed & Unassayed Multiqual®  
Liquichek™ Ethanol/Ammonia Control  
Liquichek™ Lipids Control  
Liquichek™ Microalbumin Control  
Liquichek™ Pediatric Control  
Lyphocheck® Quantitative Urine Control  
Liquichek™ & Lyphocheck® Unassayed Chemistry Controls (Human)  
Liquichek™ Urine Chemistry Control

#### Coagulation

Lyphocheck® Coagulation Control  
Liquichek™ D-dimer Control  
Lyphocheck® Hemostasis Control

#### Congenital Diseases

Liquichek™ ToRCH Plus Control  
Liquichek™ ToRCH Plus IgM Control  
VIROCLEAR® MuMZ  
VIROTROL® MuMZ  
VIROCLEAR® ToRCH  
VIROTROL® ToRCH  
VIROTROL® ToRCH-M

#### Diabetes/Hemoglobin

Liquichek™ Diabetes Control  
Lyphocheck® Diabetes Control  
Lyphocheck® Hemoglobin A1C Linearity Set  
Lyphocheck® Hemoglobin A2 Control  
Meter Trax™ Control

#### Drugs of Abuse

Liquichek™ Opiate Control  
Liquichek™ Qualitative Urine Toxicology Control  
Lyphocheck® Urine Toxicology Control  
Liquichek™ Urine Toxicology Negative Control  
Liquichek™ Urine Toxicology Controls  
Levels C1, C2, C3 & C4  
Levels C2 & C3 Low Opiate  
Levels S1, S2 & S3  
Levels S1 & S2 Low Opiate  
Levels S1E & S2E  
Levels S1E & S2E Low Opiate  
Levels S1S & S2S

#### Hematology

Liquichek™ Hematology-16 Control  
Liquichek™ Hematology-16T Control  
Liquichek™ Hematology Controls  
(A), (A-1), (C), (S) & (X)  
Liquichek™ Reticulocyte Control  
Liquichek™ Reticulocyte Controls  
(A), (A-1), (S) & (X)  
Liquichek™ Sedimentation Rate Control

#### Hepatitis & Retrovirus

Assayed VIROTROL® I-C, I-E & I-F\*  
Assayed VIROTROL® II-A & II-B\*  
VIROCLEAR®

VIROTROL® I, II, III & IV  
VIROTROL® HAV IgM  
VIROTROL® HBeAg  
VIROTROL® Hbc IgM  
VIROTROL® HIV-1 Ag\*  
VIROTROL® HIV-2

#### Immunoassay

Lyphocheck® Anemia Control  
Lyphocheck® Fertility Control  
Lyphocheck® Hypertension Markers Control  
Liquichek™ & Lyphocheck® Immunoassay Plus Controls  
Lyphocheck® Maternal Serum Control  
Liquichek™ Specialty Immunoassay Control  
Liquichek™ Tumor Marker Control  
Lyphocheck® Tumor Marker Plus Control

#### Immunology/Protein

Liquichek™ Elevated CRP Control  
Liquichek™ Immunology Control  
Lyphocheck® Immunology Plus Control  
Liquichek™ Rheumatoid Factor Control  
Liquichek™ Spinal Fluid Control

#### Molecular

AmpliTral™ III  
AmpliClear™  
AmpliPROBE™ CT/GC  
AmpliTral™ CT/GC  
AmpliTral™ HPV  
ChlamydiaPROBE™  
GonoPROBE™

#### Sexually Transmitted Diseases

VIROTROL® Syphilis Total  
VIROTROL® RPR Panel

#### Specialty

Lyphocheck® Benzo/TCA Control-Set A  
Lyphocheck® Benzo/TCA Control-Set B  
Lyphocheck® Drug Free Serum  
Lyphocheck® Endocrine Control  
Liquichek™ Serum Volatiles Control  
Lyphocheck® Urine Bone Markers Control  
Lyphocheck® Urine Metals Control  
Lyphocheck® Whole Blood Control  
Lyphocheck® Whole Blood Metals Control  
Liquichek™ Whole Blood Volatiles Control

#### Specialty Infectious Disease

CryptoTral™  
CryptoTral™ LX  
PneumoTral™  
PyloriTral™  
VIROCLEAR® EBV  
VIROTROL® Chagas  
VIROTROL® EBV  
VIROTROL® Lyme  
VIROTROL® WNV

#### Therapeutic Drug Monitoring

Liquichek™ & Lyphocheck® Therapeutic Drug Monitoring Controls (TDM)  
Lyphocheck® Whole Blood  
Immunosuppressant Control

#### Urinalysis

qUAntify® Control  
qUAntify® Plus Control  
Liquichek™ Urinalysis Control

#### Data Management Solutions

Unity™ Interlaboratory Program  
UnityConnect™  
Unity Real Time®  
Unity Real Time® online  
Unity Desktop®  
UnityWeb®  
QCNet™ (www.QCNet.com)  
WebConnect™  
Westgard Advisor™

#### External Quality Assurance

Services (EQAS)  
Clinical Chemistry Program  
Clinical Chemistry (Monthly) Program  
Hematology Program  
Hemoglobin Program  
Immunoassay Programs 1–4  
Immunoassay (Monthly) Program  
Therapeutic Drug Monitoring Program  
Urine Chemistry Program

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Axis Shield plc, Dundee, United Kingdom  
Chromogenix Instrumentation Laboratory S.p.A., Milano, Italy  
Diagnostica Stago, Inc., Parsipanny, New Jersey  
Diagon Ltd., Budapest, Hungary  
Fortress Diagnostics Ltd., Antrim, United Kingdom  
Instrumentation Laboratory Company, Lexington, Massachusetts

Roche Diagnostics, GmbH, Penzberg, Germany  
Roche Diagnostics, Indianapolis, Indiana  
Siemes Healthcare Diagnostics Inc., Newark, Delaware  
Technoclone GmbH, Vienna, Austria  
Trinity Biotech plc, Wicklow, Ireland

Bio-Rad  
Laboratories

Clinical  
Diagnostics Group

9500 Jeronimo Road  
Irvine, California 92618  
(800) 854-6737  
FAX (949) 598-1550  
bio-rad.com/qualitycontrol

Technical Service:  
(800) 854-6737

**Australia**, Bio-Rad Laboratories Pty. Ltd., Level 5, 446 Victoria Road, Gladesville NSW 2111 • Phone 61-2-9914-2800 • Telefax 61-2-9914-2888  
**Austria**, Bio-Rad Laboratories Ges.m.b.H., Hummelgasse 88/3-6, A-1130 Vienna • Phone 43-1-877-8901 • Telefax 43-1-876-5629  
**Belgium**, Bio-Rad S.A.-N.V. Begoniastraat 5, B-9810 Nazareth Eke • Phone 32-9-385-5511 • Telefax 32-9-385-6554  
**Brazil**, Bio-Rad do Brasil, Praia de Botafogo, 440-3rd Floor, Botafogo, RJ CEP 22250-040, Rio de Janeiro • Phone 5521-3237-9400 • Telefax 5521-2527-3099  
**Canada**, Bio-Rad Laboratories, Ltd., 2403 Guénette Street, Montréal, Québec H4R 2E9 • Phone 1-514-334-4372 • Telefax 1-514-334-4415  
**China**, Bio-Rad Laboratories Shanghai Ltd., 3rd Floor, #18 Dong Fang Road, Bldg E, Poly Plaza, Pudong, Shanghai, PRC 200120 • Phone 86-21-64260808 • Telefax 86-21-64264988  
**Czech Republic**, Bio-Rad spol. s r.o., Nad ostrovem 1119/7, 147 00 Prague 4 • Phone 420-241-430-532 • Telefax 420-241-431-642  
**Denmark**, Bio-Rad Laboratories, Symbion Science Park, Fruebjergvej 3, DK-2100 Copenhagen East • Phone +45-4452-1000 • Telefax +45-4452-1001  
**Finland**, Bio-Rad Laboratories, Linnanherrankuja 16, FIN-00950 Helsinki • Phone 358-9-804-22-00 • Telefax 358-9-7597-5010  
**France**, Bio-Rad, 3 boulevard Raymond Poincaré, 92430 Marnes-la-Coquette • Phone 33-1-47-95-60-00 • Telefax 33-1-47-41-91-33  
**Germany**, Bio-Rad Laboratories GmbH, Heidemannstrasse 164, D-80939 Munich • Phone +49-(0)89-318-840 • Telefax +49-(0)89-318-84100  
**Greece**, Bio-Rad Laboratories M.E.P.E., 2-4 Mesogeion Street, Fourth Floor 115 27 Athens • Phone 30-210-7774396 • Telefax 30-210-7774376  
**Hong Kong**, Bio-Rad Pacific Ltd., Unit 1101, 11/F DCH Commercial Centre, 25 Westlands Road, Quarry Bay • Phone 852-2789-3300 • Telefax 852-2789-1290  
**Hungary**, Bio-Rad Hungary Ltd., H-1082 Budapest, Futo street 47-53, Hungary • Phone +36-1-459-6100 • Telefax +36-1-459-6101  
**India**, Bio-Rad Laboratories (India) Pvt. Ltd., Bio-Rad House, 86-87, Udyog Vihar Phase IV, Gurgaon, Haryana 122 015 • Phone 1-800-180-1224 • Telefax 91-124-2398115  
**Israel**, Bio-Rad Laboratories Ltd., 14 Homa Street, New Industrial Area, Rishon Le Zion 75655 • Phone 972-3-9636050 • Telefax 972-3-9514129  
**Italy**, Bio-Rad Laboratories S.r.l., Via Cellini 18/A, 20090 Segrate, Milan • Phone +39-02-216091 • Telefax +39-02-21609-398  
**Japan**, Bio-Rad Laboratories K.K., Tennoz Central Tower 20F, 2-2-24 Higashi-Shinagawa, Shinagawa-ku, Tokyo 140-0002 • Phone 81-3-6361-7070 • Telefax 81-3-5463-8483  
**Korea**, Bio-Rad Korea Ltd., 10th Floor, HyunJuk Building, 832-41, Gangnam-gu, Seoul 135-080 • Phone 82-2-3473-4460 • Telefax 82-2-3472-7003  
**Mexico**, Bio-Rad, S.A., Avenida Eugenia 197, Piso 10-A, Col. Narvarte, C.P. 03020 Mexico, D.F. • Phone +52 (55) 5488-7670 • Telefax 52 (55) 1107-7246  
**The Netherlands**, Bio-Rad Laboratories B.V., Fokkerstraat 2-8, 3905 KV Veenendaal • Phone +31-318-540666 • Telefax +31-318-542216  
**New Zealand**, Bio-Rad New Zealand, 189 Bush Road Unit B, Albany, Auckland • Phone 64-9-415-2280 • Telefax 64-9-415-2284  
**Norway**, Bio-Rad Laboratories, Johan Scharffenbergs vei 91, N-0694 Oslo • Phone 47-23-38-41-30 • Telefax 47-23-38-41-39  
**Poland**, Bio-Rad Polska Sp. z o.o., Nakielska Str. 3, 01-106 Warsaw • Phone 48-22-3319999 • Telefax 48-22-3319988  
**Portugal**, Bio-Rad Laboratories, Lda., Edifício Prime, Ave. Quinta Grande, 53 – Fracção 3B Alfragide 26114-521 Amadora • Phone 351-21-472-7700 • Telefax 351-21-472-7777  
**Russia**, Bio-Rad Laboratorii, Business Centre "West Bridge", Leningradsky pr-t H.37A Bld. 14, 125167 Moscow • Phone 7-495-721-14-04 • Telefax 7-495-721-14-12  
**Singapore**, Bio-Rad Laboratories (Singapore) Pte. Ltd., 27 International Business Park, #01-02 iQuest @IBP, Singapore 609924 • Phone 65-6415-3170 • Telefax 65-6415-3189  
**South Africa**, Bio-Rad Laboratories (Pty) Ltd., 34 Bolton Road, Parkwood, Johannesburg 2193 • Phone 27-11-442-85-08 • Telefax 27-11-442-85-25  
**Spain**, Bio-Rad Laboratories, S.A., C/ Caléndula, 95, Edificio M. Miniparc II, El Soto de la Moraleja, 28109 Madrid • Phone 34-91-590-5200 • Telefax 34-91-590-5211  
**Sweden**, Bio-Rad Laboratories A.B., Vintergatan 1, Box 1097, S-172 22 Sundbyberg • Phone 46-8-555-127-00 • Telefax 46-8-555-127-80  
**Switzerland**, Bio-Rad Laboratories AG, Nenzlingerweg 2, CH-4153 Reinach BL • Phone 41-61-717-95-55 • Telefax 41-61-717-95-50  
**Thailand**, Bio-Rad Laboratories Ltd., 1st & 2nd Floor, Lumpini I Bldg., 239/2 Rajdamri Rd., Lumpini, Pathumwan, Bangkok 10330 • Phone 662-651-8311 • Telefax 662-651-8312  
**United Kingdom**, Bio-Rad Laboratories Ltd., Bio-Rad House, Maxted Road, Hemel Hempstead, Herts HP2 7DX • Phone +44-(0)20-8328-2000 • Telefax +44-(0)20-8328-2550

**United States**, Bio-Rad Laboratories, Irvine, CA

ECREP

**France**, Bio-Rad  
3 boulevard Raymond Poincaré  
92430 Marnes-la-Coquette  
Phone: (33) 1-4795-6000 / Fax: (33) 1-4741-9133

BIO-RAD

745

Lyphochek®

## Coagulation Control

2

A human plasma coagulation control.  
Humanplasma-Gerinnungskontrolle.  
Contrôle de coagulation du plasma humain.  
Plasma di controllo umano per coagulazione.  
Control de plasma humano para coagulación.  
Controlo de plasma humano para coagulação.  
En koagulationskontroll av humanplasma.  
En koagulationskontroll af human plasma.

Level 2

12 x 1 mL

8°C

2°C

78202

EXP 2013-06-30

UNITED STATES, Bio-Rad Laboratories, Irvine, CA  
FRANCE, Bio-Rad, Marnes-la-Coquette

5168

\*+M235Q7452/\$\$206301378202B\*

BIO-RAD

744

Lyphochek®

## Coagulation Control

1

A human plasma coagulation control.  
Humanplasma-Gerinnungskontrolle.  
Contrôle de coagulation du plasma humain.  
Plasma di controllo umano per coagulazione.  
Control de plasma humano para coagulación.  
Controlo de plasma humano para coagulação.  
En koagulationskontroll av humanplasma.  
En koagulationskontroll af human plasma.

Level 1

12 x 1 mL

8°C

2°C

78201

EXP 2013-06-30

UNITED STATES, Bio-Rad Laboratories, Irvine, CA  
FRANCE, Bio-Rad, Marnes-la-Coquette

5168

\*+M235Q7442/\$\$2063013782019\*

BIO-RAD

745X

Lyphochek®

## Coagulation Control

Trilevel  
MiniPak

3 x 1 mL

(1 per level)

A human plasma coagulation control.  
Humanplasma-Gerinnungskontrolle.  
Contrôle de coagulation du plasma humain.  
Plasma di controllo umano per coagulazione.  
Control de plasma humano para coagulación.  
Controlo de plasma humano para coagulação.  
En koagulationskontroll av humanplasma.  
En koagulationskontroll af human plasma.

8°C

2°C

78200

EXP 2013-06-30

UNITED STATES, Bio-Rad Laboratories, Irvine, CA  
FRANCE, Bio-Rad, Marnes-la-Coquette

5168

\*+M235Q745X2/\$\$2063013782008\*

BIO-RAD

746

Lyphochek®

## Coagulation Control

3

A human plasma coagulation control.  
Humanplasma-Gerinnungskontrolle.  
Contrôle de coagulation du plasma humain.  
Plasma di controllo umano per coagulazione.  
Control de plasma humano para coagulación.  
Controlo de plasma humano para coagulação.  
En koagulationskontroll av humanplasma.  
En koagulationskontroll af human plasma.

8°C

2°C

78203

EXP 2013-06-30

UNITED STATES, Bio-Rad Laboratories, Irvine, CA  
FRANCE, Bio-Rad, Marnes-la-Coquette

5168

\*+M235Q7462/\$\$206301378203D\*