Liquichek™ Blood Gas Control Levels 1, 2 and 3

REF	571	Level 1	30 x 1.7 mL
	572	Level 2	30 x 1.7 mL
	573	Level 3	30 x 1.7 mL
	570X	Trilevel MiniPak	9 x 1.7 mL

EXP 2015-03-31

21220

Level 1	21221
Level 2	21222
Level 3	21223

ENGLISH

INTENDED USE

Liquichek Blood Gas Control is an assayed quality control intended for use in monitoring the precision of an individual laboratory's measurement of pH, pCO₂ and pO₂ by blood gas instrumentation.

SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. Three levels of control are available to allow performance monitoring within the clinical range.

REAGENT

This product is a buffered bicarbonate solution in equilibrium with predetermined levels of oxygen, carbon dioxide and nitrogen.

Material Safety Data Sheet (MSDS) available for professional users on www.bio-rad.com.

STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened at 18 to 25°C. Once the control is opened, all analytes should be assayed immediately; discard the remaining material.

Avoid storage above 30°C. Freezing this control will not adversely affect performance.

This product is shipped under ambient conditions.

PROCEDURE

This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit, or reagent being used.

For optimum performance, allow a minimum of 4 hours for the ampules to equilibrate to the temperature of 20 to 23°C. Before sampling, hold the ampule between the thumb and forefinger and shake the control for 10 seconds to completely mix the contents. Restore the control to the bottom of the ampule by gentle tapping. Protect fingers with a tissue or gloves to avoid injury. Snap the ampule open and immediately introduce the control into the instrument using one of the following techniques: direct aspiration, syringe transfer or capillary mode transfer.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.

LIMITATIONS

1. This product should not be used past the expiration date.
2. This product is not intended for use as a standard.
3. This product is sensitive to most instrument related factors which cause unexpected analytical deviations. However, this control does not contain red blood cells and may not detect conditions that cause lysis of red blood cells, resulting in erroneous recovery of blood gas values.

ASSIGNMENT OF VALUES

The means and acceptable ranges reported in this insert were derived from replicate analyses and are specific for this lot of product. These tests were performed on several blood gas analyzers. For some instruments the assigned values are based on comparative studies for each make and model. Individual laboratory means should fall within the corresponding acceptable range, however, other factors such as procedural or reagent modifications may alter the means. Each laboratory should establish its own means and acceptable ranges and use the values provided only as guides.

SPECIFIC PERFORMANCE CHARACTERISTICS

This product is a stabilized liquid product manufactured under rigid quality control standards. To obtain consistent ampule-to-ampule assay values, the control requires proper storage and handling as described.

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Instrumentation Laboratory (IL), Lexington, Massachusetts
L. Eschweiler and Company, Kiel, Germany
Nova Biomedical, Waltham, Massachusetts

Radiometer of America, Westlake, Ohio
Roche Diagnostics Corporation, Indianapolis, Indiana
Siemens Healthcare Diagnostics Inc., Tarrytown, New York

Bio-Rad Laboratories comprehensive line of quality controls and QC data management solutions.

Autoimmune

Liquichek™ ANA Controls
Centromere Pattern
Homogeneous Pattern
Mitotic Spindle Pattern
Nucleolar Pattern
Speckled Pattern
Liquichek™ ANA Controls Set, Positive
(Homogeneous, Speckled, Centromere
& Nucleolar Patterns)
Liquichek™ Anti-Mitochondrial Control
Liquichek™ Anti-nDNA Control
Liquichek™ Anti-RNP Control
Liquichek™ Anti-Scl-70 Control
Liquichek™ Anti-Sm Control
Liquichek™ Anti-Smooth Muscle Control
Liquichek™ Anti-SS-A Control
Liquichek™ Anti-SS-B Control
Liquichek™ Autoimmune Negative Control

Blood Gas

Liquichek™ Blood Gas Control
Liquichek™ Blood Gas Plus CO-Oximeter
Control (Bayer 800 Series)
Liquichek™ Blood Gas Plus E Control
Liquichek™ Blood Gas Plus EGL Control

Cardiac Assessment

Liquichek™ Cardiac Markers Control LT
Liquichek™ Cardiac Markers Control LT,
Level Low
Liquichek™ Cardiac Markers Plus Control
Liquichek™ Cardiac Markers Plus Control LT
Liquichek™ CK/LD Isoenzyme Control
Liquichek™ Homocysteine Control

Chemistry

Liquichek™ Assayed Chemistry Control
Liquid Assayed & Unassayed Multiquat®
Liquichek™ Ethanol/Ammonia Control
Liquichek™ Lipids Control
Liquichek™ Microalbumin Control
Liquichek™ Pediatric Control
Liquichek™ Quantitative Urine Control
Liquichek™ & Liquichek™ Unassayed
Chemistry Controls (Human)
Liquichek™ Urine Chemistry Control

Coagulation

Liquichek™ Coagulation Control
Liquichek™ D-dimer Control
Liquichek™ Hemostasis Control

Congenital Diseases

Liquichek™ ToRCH Plus Control
Liquichek™ ToRCH Plus IgM Control
VIROCLEAR® MuMZ
VIROTROL® MuMZ
VIROCLEAR® ToRCH
VIROTROL® ToRCH
VIROTROL® ToRCH-M

Diabetes/Hemoglobin

Liquichek™ Diabetes Control
Liquichek™ Diabetes Control
Liquichek™ Hemoglobin A1C Linearity Set
Liquichek™ Hemoglobin A2 Control
Meter Trax™ Control

Drugs of Abuse

Liquichek™ Opiate Control
Liquichek™ Qualitative Urine Toxicology Control
Liquichek™ Urine Toxicology Control
Liquichek™ Urine Toxicology Negative Control
Liquichek™ Urine Toxicology Controls
Levels C1, C2, C3 & C4
Levels C2 & C3 Low Opiate
Levels S1, S2 & S3
Levels S1 & S2 Low Opiate
Levels S1E & S2E
Levels S1E & S2E Low Opiate
Levels S1S & S2S

Hematology

Liquichek™ Hematology-16 Control
Liquichek™ Hematology-16T Control
Liquichek™ Hematology Controls
(A), (A-1), (C), (S) & (X)
Liquichek™ Reticulocyte Control
Liquichek™ Reticulocyte Controls
(A), (A-1), (S) & (X)
Liquichek™ Sedimentation Rate Control

Hepatitis & Retrovirus

Assayed VIROTROL® I-C, I-E & I-F*
Assayed VIROTROL® II-A & II-B*
VIROCLEAR®

VIROTROL® I, II, III & IV
VIROTROL® HAV IgM
VIROTROL® HBsAg
VIROTROL® Hbc IgM
VIROTROL® HIV-1 Ag
VIROTROL® HIV-2

Immunoassay

Liquichek™ Anemia Control
Liquichek™ Fertility Control
Liquichek™ Hypertension Markers Control
Liquichek™ & Liquichek™ Immunoassay
Plus Controls
Liquichek™ Maternal Serum Control
Liquichek™ Specialty Immunoassay Control
Liquichek™ Tumor Marker Control
Liquichek™ Tumor Marker Plus Control

Immunology/Protein

Liquichek™ Elevated CRP Control
Liquichek™ Immunology Control
Liquichek™ Immunology Plus Control
Liquichek™ Rheumatoid Factor Control
Liquichek™ Spinal Fluid Control

Molecular

AmpliTrol™ III
AmpliClear™
AmpliPROBE™ CT/GC
AmpliTrol™ CT/GC
AmpliTrol™ HPV
ChlamydiaPROBE™
GonoPROBE™

Sexually Transmitted Diseases

VIROTROL® Syphilis Total
VIROTROL® RPR Panel

Specialty

Liquichek™ Benzo/TCA Control-Set A
Liquichek™ Benzo/TCA Control-Set B
Liquichek™ Drug Free Serum
Liquichek™ Endocrine Control
Liquichek™ Serum Volatiles Control
Liquichek™ Urine Bone Markers Control
Liquichek™ Urine Metals Control
Liquichek™ Whole Blood Control
Liquichek™ Whole Blood Metals Control
Liquichek™ Whole Blood Volatiles Control

Specialty Infectious Disease

CryptoTrol™
CryptoTrol™ LX
PneumoTrol™
PyloriTrol™
VIROCLEAR® EBV
VIROTROL® Chagas
VIROTROL® EBV
VIROTROL® Lyme
VIROTROL® WNV

Therapeutic Drug Monitoring

Liquichek™ & Liquichek™ Therapeutic
Drug Monitoring Controls (TDM)
Liquichek™ Whole Blood
Immunosuppressant Control

Urinalysis

qUAntify® Control
qUAntify® Plus Control
Liquichek™ Urinalysis Control

Data Management Solutions

Unity™ Interlaboratory Program
UnityConnect™
Unity Real Time®
Unity Real Time® online
Unity Desktop®
UnityWeb®
QCNet™ (www.QCNet.com)
WebConnect™
Westgard Advisor™

External Quality Assurance Services (EQAS)

Clinical Chemistry Program
Clinical Chemistry (Monthly) Program
Hematology Program
Hemoglobin Program
Immunoassay Programs 1-4
Immunoassay (Monthly) Program
Therapeutic Drug Monitoring Program
Urine Chemistry Program

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